

The Tao Of X-Ray Analytical Systems




Kamen Kamenov
kkamenov@issp.bas.bg



**GEORGI NADJAKOV INSTITUTE
OF SOLID STATE PHYSICS**

BULGARIAN ACADEMY OF SCIENCES



One (and only one) slide about philosophy


- Lao Tzu (Old Master) wrote “Tao Te Ching” (“The Book of the Tao and Its Power”)
- The masculine principle (Yang) - associated with strength, activity, light, hardness, and the energy of the day.
- The feminine principle (Yin) - associated with softness, receptiveness, darkness, gentleness, and the energy of the night.
- These two opposites are interdependent and complementary. Together, they represent the dynamic balance and harmony – Tao.






	XRD	XRF
Coherent scattering of X-rays	Useful signal	Noise
X-Ray florescence	Noise	Useful signal
Quality results	Crystal Lattice (Distances between atom planes)	Type of atoms
Quantitative results	Compound concentrations	Atoms concentrations
Type of measurement	Absolute	Comparative

XRD/XRF tubes



XRD Side windows


- Point focus (A line from orthogonal window)
- Anode on GND
- Cathode with -kV
- Water cooling (tape water)




XRF End window tube

- Large surface focus (30mm)
- Anode with +kV
- Cathode on GND
- Water cooling (deionized water)

X-ray tubes spectra



at $In\lambda$



XRD: line spectra, $K_{\alpha,\beta}$ shells recombination.

HV settings policy:

$$V = 5 \cdot E_{K\alpha}$$

$$I = I_{\max}$$

Example for CuK_{α}

$$V = 5 \cdot 8 \text{ kV} = 40 \text{ kV}$$

$$40 \text{ mA} < I < 55 \text{ mA}$$


70% to save tubes

XRF: continuous spectra
(**Bremsstrahlung**), accelerated electrons in anode material.

HV settings policy:

$V = V_{\max}$ or lower to excite light elements more efficiently

$$I = I_{\max} \text{ or lower}$$



Xrays - Particles or Waves



Planck Formula

$$E \cdot \lambda = h \cdot c$$

E - energy
 λ - wavelength
c - speed of light
h - Planck constant

$$E \lambda = 1.2 \text{ [kV.nm]}$$

Or


$$E \lambda = 12 \text{ [kV.Å]}$$

X-ray tube: characteristic spectra ($K_{\alpha, \beta}$) - particles
Continues spectra - waves

Diffraction: waves

Wave Dispersive X-ray
Fluorescence - WDXRF


X-Ray Diffraction - XRD




Registration: particles

Energy Dispersive X-ray
Fluorescence - EDXRF

But photon counting
detectors in XRD and WDXRD too



Bragg's Law



Bragg's Law:


Conditions to have maximum:

$$n\lambda = 2d \sin\theta$$

XRD: known λ , $d=?$

WDXRF: known d , $\lambda=?$
(but a monocrystal have to be used)

Goniometers




XRD Bragg-Bretano Geometry


- Point(line) focus x-ray source
- Divergent incident beam
- 3 Slits
- Diffraction from powder sample
- Receiving(detector) slit = slower element in the system.
- Distance Source to Sample = Sample Receiving Slit
- X-Ray Source angle = Detector angle = θ

WD-XRF Parallel Geometry

- Large spot (to 34mm) x-ray source
- Parallel incidence beam
- 2 sollers/collimators
- Diffraction from an analyzer monocrystal
- Distances are differed (detectors with collimators)
- But X-Ray Source angle = Detector angle = θ




Modern Detectors



XRD Silicon Strip Detector(SSD)

- A position sensitive solid state detector
- Every strip is equivalent of a detector slit
- $0.04^\circ - 0.01^\circ$ strip dimensions in deg.
- Efficiency $> 90\%$
- Resolution 20% to 8% ($Cu_{K\beta} - Cu_{K\alpha} / Cu_{K\alpha} = 11\%$)
Without $Cu_{K\beta}$ filter 2*


EDXRF Silicon Drift Detectors (SDD)



- Electrons drift to the central area due to the field strips.(with red)
- An integrated FET in the center with very low C
- Works from NaK_{α}
- built-in electric cooling
- 130-140 eV @ 5.9 keV MnK_{α} -Ext. 2%


$$C1U1 = C2dU2, \text{ if } C1 \ll C2, U1 \gg U2$$

Detector Electronics




XRD, WDXRF

1. Detector
2. Pulse Shape Preamplifier
3. Discriminator
4. Counter



EDXRF

1. Detector
2. Pulse Shape Preamplifier
3. Multichannel Analyzer



Critical Sample Thickness



Klenk, M.; Schenker, O.; Probst, U.; Bucher, E. X-ray fluorescence measurements of thin film chalcopyrite solar cells. Sol. Energy Mater. Sol. Cells 1999

$$I_r = C_r * K_r / \mu_r (1 - \exp(-\mu_r * \rho * d)) \quad (1)$$

When $d \rightarrow \infty$ then $I_r = C_r * K_r / \mu_r \quad (2)$

When $d \rightarrow 0$ then $I_r = C_r * K_r * \rho * d$
 $\rho * d * S = m_r \quad \text{or}$

$$I_r = m_r / S * K_r \quad (3)$$

A sample is considered infinite when its thickness is much greater than the critical depth – the layer that generates 99% of the XRF signal.

XRF is a comparative method – it works best when the sample closely matches the chemistry of the references.




For broader comparisons, mathematical alpha corrections are used to model matrix effects.

Alternatively, fused beads turn the sample into an effectively infinite thin layer, minimizing these effects from the start..

Where:

I_r [CPS]	r-element intensity
C_r [%]	r-element concentration in the sample
K_r	Constant depending of various properties of the spectrometer, tube, HV settings and depends on other elements in the sample.
μ_r [cm ² /g]	The mass absorption coefficient depends on other element in the sample.
ρ [g/cm ³]	Density of the sample
d [cm]	Sample thickness
m_r [g]	Mass of element r
S [cm ²]	Irradiated surface

TRXRF



For element r , when $d \rightarrow 0$

$$K_r = I_{0\lambda}(V, I) * P_{\lambda r} - \text{Monochromator}$$

$$I_r = m_r / S * K_r = m_r / S * I_0 * P_{\lambda r}$$

I_r [CPS] Fluorescent yield

I_0 [CPS] Initial intensity

m_r / S [g/mm²] Irradiated surface density

$P_{\lambda r}$ [mm²/g] or $\mu_{\lambda r}$ [cm²/g] (Tab. in NIST)

1. Thin film = no absorption in other elements in the sample (simple and linear calculations)


2. Monochromator = $\mu_{\lambda r}$ [cm²/g] (Tab. in NIST)

3. Grazing incidence to the substrate = Total reflection = Low background + No Air absorption = LLD ppm to θ . ppm

An Italian Story...

Carminati Marco , Borghi Giacomo.. ,32-Channel silicon strip detection module for combined X-ray fluorescence spectroscopy and X-ray diffractometry analysis, Frontiers in Physics 2022 vol10

DOI=10.3389/fphy.2022.910089




32-Channel silicon strip detection module for combined X-ray fluorescence spectroscopy and X-ray diffractometry analysis

Marco Carminati^{1,2*}, Giacomo Borghi³, Evgeny Demenev³,
Matteo Gugiatti^{1,2}, Giancarlo Pepponi³, Michele Crivellari³,
Francesco Ficorella³, Sabina Ronchin³, Nicola Zorzi³,
Evgeny Borovin⁴, Luca Lutterotti⁴ and Carlo Fiorini^{1,2}

¹DEIB, Politecnico di Milano, Milano, Italy, ²INFN Sezione di Milano, Milano, Italy, ³Fondazione Bruno Kessler (FBK), Trento, Italy, ⁴Università di Trento, Trento, Italy

- 32 silicon strip detector,
- 32 channel CUBE low noise charge integrator.
- Thermoelectric cooling module -10 C°
- 32 channel 16bit DAC with 62.5 MSps
- Better than 200 eV at 5.9 keV
- A goniometer with a radius of 15 cm
- Moderate vacuum (10 mbar)
- Mo tube operated at 50 W power

Color Diffraction
Or
Angularly Distributed Fluorescence



CaCO₃ + FeSO₄

<https://plotdigitizer.com/app>

MoK α = 17.48keV MoK β =19.61keV

Continuous spectra diffraction

6.56 - FeK α = 6,40keV +0.16

7.24 - FeK β = 7,08keV +0.18

8.23 - CuK α = 8,05keV +0.18

11.81 - BrK α =11.92keV -0.11

< 4.5 kV

CaK α = 3.69keV CaK β = 4.01keV

SK α = 2.32keV SK β = 2.46keV

Conclusions



XRD and XRF are fundamentally different techniques – one sees structure, the other sees composition – yet they share common roots in X-ray physics.

Still, human imagination knows no limits – and hybrid systems like TRXRF prove that diffraction and fluorescence can work together, delivering absolute results from a spectrometer built like a XRD machine.

Looking ahead, I believe that within the next decade, we may see the rise of color diffraction systems – a breakthrough that could redefine the entire X-ray analytics landscape.


Thanks for attention!
The End

Detectors Table




	Resolution	Efficiency	• keV Range
<u>Scintillation Detector</u>	51%-3.0keV@5.9kV	95%	>CuK α
<u>Gas Filled Proportional Counter</u>	17%-1.0keV@5.9kV	30-35%	BK α to FeK α
<u>XRD Silicon Strip Detector(SSD)</u>	20% to 8% @5.9kV	>90%	>FeKα
<u>EDXRF Silicon Drift Detectors (SDD)</u>	2%-130eV@5.9keV	>90%	>NaK α
WDXRF	0.08%- 5eV@5.9keV	< 5%	BK α


WDXRF Goniometer



Vacuum
Constant Temperature



Classic Detectors



Gas(Ar) Filled Proportional Counter

0.1 – 15keV (B, Na, Mg, Al, Si, K, Ca and Fe)


An electron-ion pair, approx. 0.03 keV

Boron (0.185 keV) produces approx. 6 pairs

30-40% efficiency

Res = 1.0keV(17%)@5.9kV

U = 1400-1900V



Scintillation Detector

8keV - ... (Cu ...)

NaI(Tl) 95% efficiency


Res = 3.0keV(51%)@5.9kV

U = 450-700V


Detector Electronics



[Figure 1-4] Circuit example of charge amplifier




KSPDC0111EA




XRD, WDXRF

1. Detector
2. Pulse Shape Preamplifier
3. Discriminator
4. Counter

[Figure 1-4] Circuit example of charge amplifier



KSPDC0111EA



EDXRF

1. Detector
2. Pulse Shape Preamplifier
(Usually integrated)
3. Multichannel Analyzer